# User Manual

***i700 wireless***  
Revision 1 (January 2022)

## CONTENTS

|     |           |
|-----|-----------|
| 02  | English   |
| 20  | العربية   |
| 38  | 中文(简体)    |
| 54  | Français  |
| 72  | Deutsch   |
| 90  | Italiano  |
| 108 | Português |
| 126 | русский   |
| 144 | Español   |
| 162 | 한국어       |

**English**

|  |    |
|--|----|
| <b>About This Guide</b> .....                              | 4  |
| <b>1 Introduction and Overview</b> .....                   | 4  |
| 1.1 Intended Use .....                                     | 4  |
| 1.2 Indication for Use .....                               | 4  |
| 1.3 Contraindications .....                                | 4  |
| 1.4 Qualifications of the Operating User .....             | 4  |
| 1.5 Symbols .....  | 4  |
| 1.6 i700 wireless Components Overview .....                | 5  |
| 1.7 Setting Up the i700 wireless System .....              | 6  |
| 1.7.1 Basic Settings of i700 wireless .....                | 6  |
| 1.7.2 How to Use the Wireless Hub .....                    | 7  |
| 1.7.3 How to Use the Battery .....                         | 7  |
| 1.7.4 How to Use the Battery Charger .....                 | 7  |
| 1.7.5 How to Mount the Handle .....                        | 8  |
| 1.7.6 How to Place on the Desktop Cradle .....             | 8  |
| 1.7.7 How to Install the Wall Mount Holder .....           | 8  |
| <b>2 Medit Scan for Clinics Overview</b> .....             | 8  |
| 2.1 Introduction .....                                     | 8  |
| 2.2 Installation .....                                     | 8  |
| 2.2.1 System Requirements .....                            | 8  |
| 2.2.2 Medit Scan for Clinics Installation Guide .....      | 9  |
| 2.2.3 Medit Scan for Clinics User Guide .....              | 9  |
| <b>3 Maintenance</b> .....                                 | 10 |
| 3.1 Calibration .....                                      | 10 |
| 3.1.1 How to Calibrate i700 wireless .....                 | 10 |
| 3.1.2 Auto Calibration Tool (Sold Separately) .....        | 10 |
| 3.2 Cleaning, Disinfection, Sterilization Procedure .....  | 10 |
| 3.2.1 Reusable Tip – Sterilization .....                   | 10 |
| 3.2.2 Reusable Tip – Cleaning and Disinfection .....       | 10 |
| 3.2.3 Mirror .....   | 10 |
| 3.2.4 Handpiece .....                                      | 11 |
| 3.2.5 Other Components .....                               | 11 |
| 3.3 Disposal .....   | 11 |
| 3.4 Battery Storage .....                                  | 11 |
| 3.5 The Battery Usage Precautions and Disposal Guide ..... | 11 |
| 3.6 Updates in Medit Scan for Clinics .....                | 12 |
| <b>4 Safety Guide</b> .....                                | 12 |
| 4.1 System Basics .....                                    | 12 |
| 4.2 Proper Training .....                                  | 12 |
| 4.3 In Case of Equipment Failure .....                     | 13 |
| 4.4 Hygiene .....  | 13 |
| 4.5 Electrical Safety .....                                | 13 |
| 4.6 Eye Safety .....                                       | 13 |
| 4.7 Explosion Hazards .....                                | 13 |
| 4.8 Pacemaker and ICD Interference Risk .....              | 13 |
| <b>5 Electro-Magnetic Compatibility Information</b> .....  | 14 |
| 5.1 Electromagnetic Emissions .....                        | 14 |
| 5.2 Electromagnetic Immunity .....                         | 14 |
| <b>6 Wireless Compliance Information</b> .....             | 16 |
| 6.1 FCC Compliance Statement .....                         | 16 |
| 6.2 IC Compliance Statement .....                          | 16 |
| 6.3 KC Compliance Statement .....                          | 17 |
| 6.4 TELEC (Japan) Compliance Statement .....               | 17 |
| <b>7 Specifications</b> .....                              | 17 |

## About this guide

### Convention in This Guide

This Guide uses various symbols to highlight important information to ensure correct usage, prevent injury to the user and others, and prevent property damage. The meanings of the symbols used are described below.

The WARNING symbol indicates information that, if ignored, could result in a medium risk of personal injury.

The CAUTION symbol indicates safety information that, if ignored, could result in a slight risk of personal injury, property damage, or damage to the system.

The TIPS symbol indicates hints, tips, and additional information for optimal operation of the system.

## 1. Introduction and Overview

### 1.1 Intended Use

The i700 wireless system is an intraoral 3D scanner intended to record topographical characteristics of teeth and surrounding tissues digitally. The i700 wireless system produces 3D scans for use in computer-assisted design and manufacturing of dental restorations.

### 1.2 Indication for Use

The i700 wireless system should be used on patients who require 3D scanning for dental treatments such as:

- Single Custom Abutment
- Inlays & Onlays
- Single Crown
- Veneer
- 3-Unit Implant Bridge
- Up to 5-Unit Bridge
- Orthodontics
- Implant Guide
- Diagnosis Model

The i700 wireless system can also be used in full-arch scans, but various factors (intraoral environment, operator's expertise, and laboratory workflow) may affect the final results.

### 1.3 Contraindications

The i700 wireless system is not intended to be used to create images of the internal structure of teeth or the supporting skeletal structure.

### 1.4 Qualifications of the Operating User

- The i700 wireless system is designed for use by individuals with professional knowledge in dentistry and dental laboratory technology.
- The user of the i700 wireless system is solely responsible for determining whether or not this device is suitable for a particular patient's case and circumstances.
- The user is solely responsible for the accuracy, completeness, and adequacy of all data entered the i700 wireless system and the provided software. The user should check the accuracy of the results and assess each individual case.
- The i700 wireless system must be used in accordance with its accompanying User Guide.
- Improper use or handling of the i700 wireless system will void its warranty. If you require additional information on the proper use of the i700 Wireless system, please contact your local distributor.
- The user is not allowed to modify the i700 wireless system.

### 1.5 Symbols

| No | Symbol  | Description   |
|----|---|---|
| 1  |  | Serial Number of the Object                         |
| 2  |  | Medical Device                                      |
| 3  |  | Date of Manufacture                                 |
| 4  |  | Manufacturer  |
| 5  |  | Caution   |
| 6  |  | Warning   |
| 7  |  | Instructions for User Manual                        |
| 8  |  | The Official Mark of Europe Certificate             |
| 9  |  | Authorized Representative in the European Community |
| 10 |  | BF Type of Applied Part                             |
| 11 |  | WEEE Mark   |
| 12 |  | Prescription Use (U.S.A.)                           |
| 13 |  | MET Mark  |
| 14 |  | AC  |
| 15 |  | DC  |

|    |   |                                 |
|----|---|---------------------------------|
| 16 |  | Temperature Limitation          |
| 17 |  | Humidity Limitation             |
| 18 |  | Atmospheric Pressure Limitation |
| 19 |  | Fragile                         |
| 20 |  | Keep Dry                        |
| 21 |  | This Way Up                     |
| 22 |  | Seven-Layer Stacking Prohibited |
| 23 |  | Consult Instructions For Use    |

## 1.6 i700 wireless Components Overview

| No | Item   | Qty                       | Appearance  |
|----|--|---------------------------|---|
| 1  | i700 wireless Handpiece                                    | 1ea                       |    |
| 2  | Wireless Hub   | 1ea                       |    |
| 3  | Rechargeable Battery                                       | 3ea                       |    |
| 4  | Extension Cable Battery                                    | 1ea                       |    |
| 5  | Battery Charger  | 1ea                       |    |
| 6  | Handpiece Cover  | 1ea                       |    |
| 7  | Reusable Tip   | 4ea                       |   |
| 8  | Handle   | 1ea                       |  |
| 9  | Calibration Tool   | 1ea                       |  |
| 10 | Practice Model   | 1ea                       |  |
| 11 | Wrist Strap  | 1ea                       |  |
| 12 | Desktop Cradle   | 1ea                       |  |
| 13 | Wall Mount Holder  | 1ea                       |  |
| 14 | USB 3.0 Cable (C to A)                                     | 1ea                       |  |
| 15 | Power Delivery Cable (C to C)                              | 1ea                       |  |
| 16 | Medical Adapter for Wireless Hub                           | 1ea                       |  |
| 17 | Medical Adapter for Battery Charger<br>(Same as above)     | 1ea                       |  |
| 18 | Power Cord   | 2ea                       |  |
| 19 | USB Flash Drive<br>(Installer of Medit Scan for Clinics)   | 1ea                       |  |
| 20 | User Guide   | 1ea                       |  |
| 21 | Auto Calibration Tool<br>(1 rechargeable battery included) | 1ea<br>(*Sold Separately) |  |

**CAUTION**

- Keep the practice model in a cool place away from direct sunlight. A discolored practice model may affect the results of the practice mode.
- The strap is designed specifically for the weight of the i700 wireless and should not be used with other products.
- Medit Scan for Clinics is included in the USB drive. This product is optimized for PC, and using other devices is not recommended. Do not use anything other than a USB port. It may cause malfunction or fire.
- The auto-calibration tool is not included in the i700 wireless package. It can be bought separately from your place of purchase.

**1.7 Setting Up the i700 wireless System**

**1.7.1 Basic Settings of i700 wireless**

**Connect the i700 wireless Wireless Hub (1)**

- ① Connect the USB 3.0 cable (C to A) to the wireless hub.      ② Connect the adapter to the wireless hub.

- ③ Connecting the power cord to the adapter automatically powers the wireless hub.

- ④ Connect the USB 3.0 cable's A port connector to a PC.

※ The USB port is used only for signal transmission.

**Connect the i700 wireless Wireless Hub (2)**

- ① Connect the power delivery cable (C to C) to the wireless hub.      ② Connect the power delivery cable to a PC.

※ 9 Vdc power is supplied to the USB port.

**Turning On the i700 wireless**

- ① Insert the battery into the i700 wireless handpiece and press the power button.      ② When the power is supplied, the LED on the upper part of the i700 wireless handpiece turns blue.

- ③ Three LEDs on the bottom of the i700 wireless handpiece indicate the remaining battery level.

- 3 lights : 80 - 100%
- 2 lights : 50 - 79%
- 1 light : 20 - 49%
- 1 blinking light : 1 - 19%
- LED lights turned off : 0%

#### Turning Off the i700 wireless

Press and hold the power button at the bottom of the i700 wireless handpiece for 3 seconds. If the rechargeable battery is removed without turning the device off, it may shorten the lifespan of the i700 wireless and the battery.

#### 1.7.2 How to Use the Wireless Hub

The i700 wireless is a wireless device that operates through a wireless module. For this reason, the i700 wireless handpiece has a transmitter and the wireless hub has a receiver. The i700 wireless system uses two types of frequencies to transmit data and control the i700 wireless handpiece. Power is supplied by connecting the adapter cable to the wireless hub power port. The power turns off when the adapter port is removed.

When the i700 wireless is turned on, it will automatically try to pair with the wireless hub. Both the i700 wireless handpiece and wireless hub must be powered and placed in a range of each other to pair. When pairing is in progress, the LED on top of the i700 wireless handpiece flashes. When pairing is complete, the LED lights up. You can start scanning when the pairing is complete.

- i700 wireless uses two antenna modules : 60 GHz for data transmission and 2.4 GHz for the control.
- The actual frequency depends on the local regulations.
- The actual operational range is up to 5 m, and it may vary depending on the environment.
- 60 GHz Frequency : 57 - 64 GHz
- 2.4 GHz Frequency : 2.4 - 2.5 GHz

Power is supplied from the connected PC without using the adapters when connected with a power delivery cable. In this case, the PC must be turned on. If the power delivery cable is removed from the PC, the wireless hub will be automatically turned off and the connection status, such as pairing, will be initialized.

#### 1.7.3 How to Use the Battery

- Rechargeable Battery
  - Li-ion
  - 3.6 V, 3100 mAh, 11.6 Wh
  - Insert the battery at the bottom of the i700 wireless handpiece. Insert the battery into the i700 wireless handpiece with the terminal aligned correctly.

- Battery life lasts up to 1 hour.  
Battery life may vary depending on the user environment and the battery cycle count.

- Extension Cable Battery
  - The extension cable battery is a cable-type wired battery that is used instead of a rechargeable battery. No charging is required as power is supplied through the cable.
  - Connect the battery-shaped connector to the i700 wireless handpiece and connect the cable to the battery charger terminal.

- Press the power button on the bottom of the i700 wireless handpiece to supply power.

#### 1.7.4 How to Use the Battery Charger

- Connect the adapter port to the battery charger power port to supply power. The power turns off when the adapter port is removed.

- Insert the rechargeable battery into the battery charger with the charging terminals oriented correctly.

- The battery charger is only for rechargeable batteries. It takes up to 2 hours and 30 minutes to be fully charged and the charging time may vary depending on the user environment and the battery cycle count.

The LED on the charger blinks in blue when the battery is charging. When fully charged, the LED lights up in blue.

If the battery is not correctly inserted into the battery charger, the LED on the charger will blink in red. In this case, remove the battery from the charger, check both terminals of the battery and the charger for foreign objects, wipe them gently with a soft cloth, and then insert the battery again.

### 1.7.5 How to Mount the Handle

The i700 wireless body is equipped with a wireless signal transmitter located where the i700 wireless logo is. Depending on your experiences and habits, you may hold the area where the transmitter is mounted. Covering the transmitter area may interfere with the wireless hub's communication. Therefore, a handle is provided to mount on the i700 wireless handpiece for a more comfortable hand-gripping method.

- Turn over the i700 wireless handpiece to find the silicone body.
- Remove the silicone with your hand.

- Fully fasten the bolts of the handle to the handle mounting hole on the i700 wireless body.
- Tighten it clockwise using the knob on the handle.

- You can now use it by holding the handle. If you want to remove the handle, proceed in the reverse order of this instruction.

### 1.7.6 How to Place on the Desktop Cradle

Without the Handle

With the Handle

### 1.7.7 How to Install the Wall Mount Holder

## 2. Medit Scan for Clinics Overview

### 2.1 Introduction

Medit Scan for Clinics provides a user-friendly working interface to digitally record topographical characteristics of teeth and surrounding tissues using the i700 wireless system.

### 2.2 Installation

#### 2.2.1 System Requirements

## Minimum System Requirements

|         | Laptop  | Desktop                |
|---------|---|------------------------|
| CPU     | Intel Core i7 - 10750H  | Intel Core i7 - 10700K |
|         | AMD Ryzen 7 4800H   | AMD Ryzen 7 3800X      |
| RAM     | 32 GB   |                        |
| Graphic | NVIDIA GeForce RTX 1660/2060/3060 (Above 6GB)<br>AMD Radeon is not supported. |                        |
| OS      | Windows 10 Pro 64-bit<br>Windows 11 Pro 64-bit                                |                        |

## Recommend System Requirements

|         | Laptop  | Desktop                |
|---------|---|------------------------|
| CPU     | Intel Core i7 - 11800H  | Intel Core i7 - 11700K |
|         | Intel Core i7 - 12700H  | Intel Core i7 - 12700K |
|         | AMD Ryzen 7 5800H   | AMD Ryzen 7 5800X      |
|         | AMD Ryzen 9 5900H   |                        |
| RAM     | 32 GB   |                        |
| Graphic | NVIDIA GeForce RTX 2070/2080/3070/3080/3090 (Above 8GB)<br>AMD Radeon is not supported. |                        |
| OS      | Windows 10 Pro 64-bit<br>Windows 11 Pro 64-bit  |                        |

For accurate and up-to-date system requirements, please visit [www.meditlink.com](http://www.meditlink.com).

Use PC and monitor certified IEC 60950, IEC 55032, IEC 55024.

The device may not work when using cables other than the USB 3.0 cable provided by Medit. Medit is not responsible for any problems caused by cables other than the USB 3.0 cable provided by Medit. Be sure to use only the USB 3.0 cable included in the package.

## 2.2.2 Medit Scan for Clinics Installation Guide


- ① Run the "Medit\_Scan\_for\_Clinics\_XXX.exe" file.

- ② Select the setup language and click "Next."


- ③ Select the installation path.

- ④ Read the "License Agreement" carefully before checking "I agree to the License terms and conditions." and then click "Install."

- ⑤ It may take several minutes to finish the installation process. Please do not shut down the PC until the installation is complete.

- ⑥ After the installation is complete, restart the PC to ensure optimal program operation.

The installation will not be processed while the i700 wireless system is connected to a PC. Please be sure to disconnect the i700 wireless USB 3.0 cable from the PC before the installation.

### 2.2.3 Medit Scan for Clinics User Guide

Please refer to the User Guide of Medit Scan for Clinics: Medit Scan for Clinics > Menu > User Guide.

#### CAUTION

- Equipment maintenance should only be carried out by a Medit employee or Medit-Certified company or personnel.
- In general, users are not required to perform maintenance work on the i700 wireless system besides calibration, cleaning, and sterilization. Preventive inspections and other regular maintenance are not required.

## 3. Maintenance

### 3.1 Calibration

Periodic calibration is required to produce precise 3D models. You should perform calibration when :

- The quality of the 3D model is not reliable or accurate when compared to previous results.
  - Environmental conditions such as temperature have changed.
  - The calibration period has expired.
- You can set the calibration period in the Menu > Settings > Calibration Period (Days).

The calibration panel is a delicate component.

Do not touch the panel directly. Check the calibration panel if the calibration process is not performed properly. If the calibration panel is contaminated, please contact your service provider.

The calibration panel is a delicate component.

You can set the calibration period via Menu > Settings > Calibration Period (Days). The default calibration period is 14 days.

#### 3.1.1 How to Calibrate i700 wireless

- ① Turn on the i700 wireless and launch the Medit Scan for Clinics.
- ② Run the Calibration Wizard from Menu > Settings > Calibration.
- ③ Prepare the calibration tool and the i700 wireless handpiece.
- ④ Turn the dial of the calibration tool to position **1**.
- ⑤ Put the i700 wireless handpiece into the calibration tool.
- ⑥ Click "Next" to start the calibration process.
- ⑦ When the calibration tool is mounted properly in the correct position **1**, the system will automatically acquire the data.
- ⑧ When data acquisition is complete at position **1**, turn the dial to the next position.
- ⑨ Repeat the steps for positions **2** – **8** and the **LAST** position.
- ⑩ When data acquisition is complete at the **LAST** position, the system will automatically calculate and show the calibration results.

#### 3.1.2 Auto Calibration Tool (Sold Separately)

The i700 wireless auto calibration tool accessory can be purchased separately. This handy auto calibration tool will automatically perform calibration by calibrating the i700 wireless handpiece without having to turn the calibration dial. Please refer to Medit Scan for Clinics for details.

### 3.2 Cleaning, Disinfection, Sterilization Procedure

#### 3.2.1 Reusable Tip – Sterilization

The reusable tip is the part that is inserted into the patient's mouth during scanning. The tip is reusable for a limited number of times. The tip needs to be cleaned and sterilized between patients to avoid cross-contamination.

- The tip should be cleaned manually using a disinfecting solution. After cleaning and disinfecting, inspect the mirror inside the tip to ensure there aren't any stains or smudges.
- Repeat the cleaning and disinfection process if necessary. Carefully dry the mirror using a paper towel.
- Insert the tip into a paper sterilization pouch and seal it, making sure that it is airtight. Use either a self-adhesive or heat-sealed pouch.
- Sterilize the wrapped tip in an autoclave with the following conditions:
  - Sterilize for 30 minutes at 121°C (249.8°F) at gravity type and dry for 15 minutes.
  - Sterilize for 10 minutes at 135°C (275°F) at gravity type and dry for 30 minutes.
  - Sterilize for 4 minutes at 134°C (273.2°F) at pre-vacuum type and dry for 20 minutes.
- Use an autoclave program that dries the wrapped tip before opening the autoclave.
- Scanner tips can be re-sterilized up to 150 times and thereafter must be disposed of as described in the disposal section.
- Autoclave times and temperatures may vary depending on the autoclave type and manufacturer. For this reason, it may not be able to meet the maximum number of times. Please refer to the user's manual of the autoclave manufacturer you are using to determine whether the required conditions are met.

#### 3.2.2 Reusable Tip – Cleaning and Disinfection

- Clean the tip immediately after use with soap water and a brush. We recommend using a mild dishwashing liquid. Make sure the mirror of the tip is completely clean and stain-free after cleaning. If the mirror appears stained or foggy, repeat the cleaning process and rinse thoroughly with water. Dry the mirror carefully with a paper towel.
- Clean with Caviwipes according to the following conditions. Please refer to the Caviwipes instruction manual for proper use.
  - Caviwipes : disinfect for 3 minutes and dry for 5 minutes.
  - Caviwipes-1 : disinfect for 1 minute and dry for 5 minutes.
  - Caviwipes-2 : disinfect for 2 minutes and dry for 5 minutes.
- Disinfect the tip using Wavicide-01 for 45 to 60 minutes. Rinse the tip thoroughly. Please refer to the Wavicide-01 Solution instruction manual for proper use.
- Remove the tip from the used solution and rinse thoroughly after cleaning and sterilization.
- Use a sterilized and non-abrasive cloth to dry the mirror and the tip gently.

**⚠ CAUTION**

- The mirror found in the tip is a delicate optical component that should be handled with care to ensure optimal scan quality. Be careful not to scratch or smudge it as any damage or blemishes may affect the data acquired.
- Make sure to always wrap the tip before autoclaving. If you autoclave an exposed tip, it will cause stains on the mirror, which cannot be removed. Check the autoclave manual for more information.
- Tips that have been cleaned, disinfected, and sterilized must remain sterile until they are used on the patient.
- Medit is not responsible for any damage, such as distortion of the tip, that occurs during cleaning, disinfection, or sterilization operations that are not following the guidelines above.

**3.2.3 Mirror**

The presence of impurities or smudges on the tip mirror may lead to poor scan quality and an overall poor scanning experience. In such situations, clean the mirror following the steps below:

- ① Disconnect the scanner tip from the i700 wireless handpiece.
- ② Pour alcohol on a clean cloth or cotton-tipped swab and wipe the mirror. Make sure to use alcohol that is free of impurities or it may stain the mirror. You can use either ethanol or propanol (ethyl-/propyl alcohol).
- ③ Wipe the mirror dry using a dry, lint-free cloth.
- ④ Make sure the mirror is free of dust and fibers. Repeat the cleaning process as necessary.

**3.2.4 Handpiece**

After treatment, clean and disinfect all other surfaces of the i700 wireless handpiece except for the scanner front (optical window) and end (air vent hole). Cleaning and disinfecting must be done with the device turned off. Use the device only after it is completely dry.

The recommended cleaning and disinfecting solution is denatured alcohol (ethyl alcohol or ethanol) – typically 60 – 70% Alc/Vol.

The general cleaning and disinfecting procedures are as follows:

- ① Turn off the device using the power button.
- ② Unplug all the cables from the wireless hub.
- ③ Clean the filter on the front end of the i700 wireless handpiece.
  - If alcohol is poured directly into the filter, it may seep inside the i700 wireless handpiece and cause a malfunction.
  - Do not clean the filter by pouring alcohol or cleaning solution directly into the filter. The filter must be gently wiped with a cotton or soft cloth moistened with alcohol. Do not wipe by hand or apply excessive force.
  - Medit is not responsible for any damage or malfunction that occurs during cleaning that does not follow the guidelines above.
- ④ Put the cover on the front of the i700 wireless handpiece after cleaning the filter.
- ⑤ Pour the disinfectant onto a soft, lint-free, and non-abrasive cloth.
- ⑥ Wipe the scanner surface with the cloth.
- ⑦ Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

**⚠ CAUTION**

- Do not clean the i700 wireless handpiece when the device is turned on as the fluid may enter the scanner and cause malfunction.
- Use the device after it is completely dry.
- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

**3.2.5 Other Components**

- Pour the cleaning and disinfecting solution onto a soft, lint-free, and non-abrasive cloth.
- Wipe the component surface with the cloth.
- Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

**⚠ CAUTION**

- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

**3.3 Disposal****⚠ CAUTION**

- The scanner tip must be sterilized before disposal. Sterilize the tip as described in the section "3.2.1 Reusable Tip – Sterilization."
- Dispose of the scanner tip as you would any other clinical waste.
- Other components are designed to conform with the following directives:
  - RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment. (2011/65/EU)
  - WEEE, Waste Electrical and Electronic Equipment Directive. (2012/19/EU)

**3.4 Battery Storage**

- Put it in a pack or a box and store indoors in a cool environment out of direct sunlight.
- Store the battery in a dry place with an ambient temperature of from -20°C to +30°C (-4°F to +86°F).
- If not used for a long time, the battery self-discharge may accelerate and slip into sleep mode. To minimize the effect of deactivation, store the battery pack in between +10°C - +30°C (+50°F to +86°F).
- When charging for the first time after long-term storage, the battery capacity may decrease due to pack deactivation. Restore the battery through several full charge/discharge cycles.
- If the battery is stored for more than 6 months, it should be charged at least once every 6 months to prevent the decrease in shelf life due to self-discharge.

**⚠ CAUTIONS: Battery Safety Specifications**

| Safety Specifications |                |                 |
|-----------------------|----------------|-----------------|
| Overcharge            | Sense Voltage  | 4.225 V ± 0.020 |
|                       | Steady Voltage | 4.025 V ± 0.03  |
|                       | Sense Delay    | 1.0 s ± 0.2     |
| Over-Discharge        | Sense Voltage  | 2.50 V ± 0.035  |
|                       | Steady Voltage | 2.90 V ± 0.50   |
|                       | Sense Delay    | 64 ms ± 12.8    |

|                                       |                   |                     |
|---------------------------------------|-------------------|---------------------|
| Overcurrent                           | Sense (Charge)    | 10.0 A +5.0 / -4.0  |
|                                       | Sense Delay       | 8.0 ms ± 1.6        |
|                                       | Sense (Discharge) | 10.0 A + 4.4 / -3.8 |
|                                       | Sense Delay       | 8.0 ms ± 1.6        |
| Current consumption at operating mode |                   | Max. 150.0 µA       |

The possible safety specifications are determined by the propulsion control module (PCM) in the parts list.

### 3.5 The Battery Usage Precautions and Disposal Guide

#### CAUTION

- Be sure to fully understand how to replace the battery before use.
  - Use a charger suitable for the specified voltage and current.
  - Do not attempt to reverse the charge. Reverse charging may increase the gas pressure inside the battery and cause the battery to leak.
  - Do not attempt to recharge a fully charged battery. Repeated overcharging could result in battery performance degradation and overheating.
  - Charging efficiency decreases at temperatures above +40°C (+104°F).
  - Do not short-circuit the positive (+) and negative (-) terminals with metallic objects such as metal wires, necklaces, or chains.
  - To avoid malfunction or damage, do not drop or throw the battery.
- 
- Do not deform a battery by excessive pressure.
  - Do not directly solder anything onto the battery pack.
  - Do not allow children to change batteries without adult supervision.
  - Do not dispose of the battery as general waste and separate it from recyclables.
  - Do not dispose of or throw the battery into a fire. Heat can cause a battery explosion and a fire.
  - Separate them from each other when disposing of the secondary batteries with different electrochemical systems.
  - Dispose of the battery by discharging it to prevent heat from short-circuiting.
  - Battery disposal methods may vary by country and region. Dispose of used batteries under local laws and regulations.

### 3.6 Updates in Medit Scan for Clinics

Medit Scan for Clinics automatically checks for updates when the software is in operation. If a new version of the software is released, the system will automatically download it.

## 4. Safety Guide

Please adhere to all the safety procedures as detailed in this User Guide to prevent human injury and equipment damage. This document uses the words WARNING and CAUTION when highlighting precautionary messages.

Carefully read and understand the guidelines, including all preventive messages as prefaced by the words WARNING and CAUTION. To avoid bodily injury or equipment damage, make sure to adhere strictly to the safety guidelines. All instructions and precautions as specified in the Safety Guide must be observed to ensure the system's proper functionality and personal safety.

The i700 wireless system should only be operated by dental professionals and technicians who are trained to use the system. Using the i700 wireless system for any purpose other than its intended usage as outlined in the section "1.1 Intended Use" may result in injury or damage to the equipment. Please handle the i700 wireless system according to the guidelines in the safety guide.

### 4.1 System Basics

The i700 wireless system is a high-precision optical measuring device. Acquaint all the following safety and operating instructions before the installation, usage, and operation of i700 wireless.

#### CAUTION

- The USB 3.0 cable connected to the wireless hub is the same as a standard USB cable connector. However, the device may not operate normally if a standard 3.0 USB cable is used with i700 wireless.
- The wireless hub is designed specifically for the i700 wireless and should not be used with any other device.
- This device has been assessed for conformity for use in business environments. When used in residential environments, there is a risk of radio wave interference.
- If the product has been stored in a cold environment, give it time to adjust to the temperature of the environment before use. If used immediately, condensation may occur, which may damage the electronic parts inside the unit.
- Ensure that all components provided are free from physical damage. Safety cannot be guaranteed if there is any physical damage to the unit.
- Before using the system, check that there are no issues such as physical damage or loose parts. If there is any visible damage, do not use the product and contact the manufacturer or your local representative.
- Check the i700 wireless handpiece and its accessories for any sharp edges.
- When not in use, the i700 wireless system should be kept mounted on a desk stand or a wall mount stand.
- Do not install the desk stand on an inclined surface.
- Do not place any object on the i700 wireless system.
- Do not place the i700 wireless system on any heated or wet surface.
- Do not block the air vents located at the rear of the i700 wireless system. If the equipment overheats, the i700 wireless system may malfunction or stop working.
- The battery of the i700 wireless can only be compatible with the i700 wireless system.
- Do not touch the charging terminals of the rechargeable battery with your hands or other instruments.
- If the charging terminal of the rechargeable battery is damaged, do not use it and contact the manufacturer or regional manager.
- If the shape of the rechargeable battery has been deformed by dropping or physical damage, never use it, and contact the manufacturer or regional manager.
- Extension cable battery is not designed for charging. Do not attempt to charge it with the battery charger.
- Use only the extension cable battery provided by the manufacturer.
- Do not spill any liquid on the i700 wireless system.
- The i700 wireless handpiece and other included components are made of electronic components. Do not allow any kind of liquid or foreign objects to enter.
- Do not pull or bend the cable connected to the i700 wireless system.
- Carefully arrange all the cables so that you or your patient do not trip or get caught in the cables. Any pulling tension on the cables may cause damage to the i700 wireless system.
- Always place the plug of the power cord of the i700 wireless system in an easily accessible location.
- Always keep an eye on the product and your patient while using the product to check for abnormalities.
- Proceed with calibration, cleaning, disinfection, and sterilization in accordance with the contents of the user guide.
- If you drop the i700 wireless tip on the floor, do not attempt to reuse it. Discard the tip immediately as there is a risk that the mirror attached to the tip may have been dislodged.
- Due to its fragile nature, the i700 wireless tips should be handled with care. To prevent damage to the tip and its internal mirror, be

careful to avoid contact with the patient's teeth or restorations.

- If the i700 wireless system is dropped on the floor or if the unit is impacted, it must be calibrated before use. If the instrument is unable to connect to the software, consult the manufacturer or authorized resellers.
- If the equipment fails to operate normally, such as having issues with accuracy, stop using the product, and contact the manufacturer or authorized resellers.
- Install and use only approved programs to ensure the proper functionality of the i700 wireless system.
- In the event of a severe accident involving the i700 wireless system, notify the manufacturer and report it to the competent national authority of the country where the user and patient reside.
- If the PC with the software installed does not have security software or if there is a risk of malicious code intrusion into the network, the PC may be breached with malware (malicious software such as viruses or worms that damage your computer).
- The software for this product must be used in compliance with medical and personal information protection laws.

#### 4.2 Proper Training

##### WARNING

Before using your i700 wireless system on patients :

- You should have been trained to use the system or read and fully understand this User Guide.
- You should be familiar with the safe use of the i700 wireless system, as detailed in this User Guide.
- Before use or after changing any settings, the user should check that the live image is displayed properly in the camera preview window of the program.

#### 4.3 In Case of Equipment Failure

##### WARNING

If your i700 wireless system is not working properly or if you suspect that there is a problem with the equipment :

- Remove the device from the patient's mouth and discontinue use immediately.
- Disconnect the device from the PC and check for errors.
- Remove the rechargeable battery from the i700 wireless system.
- Contact the manufacturer or authorized resellers.
- Modifications to the i700 wireless system are prohibited by law as they may compromise the safety of the user, patient, or a third party.

#### 4.4 Hygiene

##### WARNING

For clean working conditions and patient safety, ALWAYS wear clean surgical gloves when :

- Handling and replacing the tip.
- Using the i700 wireless system on patients.
- Touching the i700 wireless system.

##### WARNING

The i700 wireless system and its optical window should be always kept clean. Before using the i700 wireless system on a patient, be sure to :

- Sterilize the i700 wireless system as described in the section "3.2 Cleaning, Disinfection, Sterilization Procedure."
- Use a sterilized tip.

#### 4.5 Electrical Safety

##### WARNING

- The i700 wireless system is a Class I device. The i700 wireless system includes the i700 wireless handpiece, wireless hub, battery charger, and rechargeable battery altogether.
- To prevent electric shock, the i700 wireless system must only be connected to a power source with a protective earth connection. If you cannot insert the i700 wireless-supplied plug into the main outlet, contact a qualified electrician to replace the plug or outlet. Do not try to circumvent these safety guidelines.
- Do not use a grounding-type plug connected to the i700 wireless system for any other purpose than its intended use.
- The i700 wireless system only uses RF energy internally. The amount of RF radiation is low and does not interfere with surrounding electromagnetic radiation.
- There is a risk of electric shock if you attempt to access the inside of the i700 wireless system. Only qualified service personnel should access the system.
- Do not connect the i700 wireless system to a regular power strip or extension cord, as these connections are not as safe as grounded outlets. Failure to adhere to these safety guidelines may result in the following hazards :
  - All connected equipment's total short circuit current may exceed the limit specified in EN/IEC 60601-1.
  - The impedance of the ground connection may exceed the limit specified in EN/IEC 60601-1.
- Do not place liquids such as beverages near the i700 wireless system and avoid spilling any liquid on the system.
- Never spill liquid of any kind on the i700 wireless system.
- Condensation due to changes in temperature or humidity can cause moisture buildup inside the i700 Wireless system, which may damage the system. Before connecting the i700 wireless system to a power supply, be sure to keep the i700 wireless system at room temperature for at least two hours to prevent condensation. If condensation is visible on the product surface, the i700 wireless should be left at room temperature for more than 8 hours.
- You should only disconnect the i700 wireless system from the power supply via its power cord or rechargeable battery.
- When disconnecting the power cord or the rechargeable battery, hold the surface to remove it.
- Before disconnecting, make sure to turn off the power on the device using the power switch on the handpiece.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If used in a residential environment (for which CISPR 11 Class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services.
- Use only batteries provided for use with the i700 wireless. Other batteries could damage the i700 wireless system.
- Avoid pulling on the communication cables, power cables, etc. used with the i700 wireless system.
- Use only medical adapters provided for use with the i700 wireless. Other adapters could damage the i700 wireless system.
- Do not touch the connectors of the device and the patient simultaneously.

#### 4.6 Eye Safety

##### WARNING

- The i700 wireless system projects a bright light from its tip during scanning.
- The bright light projected from the tip of the i700 wireless is not harmful to the eyes. However, you should not look directly at the bright light nor aim the light beam into the eyes of others. Generally, intense light sources can cause eyes to become brittle and the likelihood of secondary exposure is high. As with other intense light source exposure, you may experience a temporary reduction in visual acuity, pain, discomfort, or visual impairment, all of which increase the risk of secondary accidents.

## 5. Electro-Magnetic Compatibility Information

- There is an LED that emits UV-C wavelengths inside the i700 wireless handpiece. It is irradiated only inside the i700 wireless handpiece and does not go outside. The blue light visible inside the i700 wireless handpiece is for guidance, not UV-C light. It is harmless to the human body.
- The UV-C LED operates with a wavelength of 270 - 285 nm.
- Disclaimer for Risks Involving Patients with Epilepsy  
The Medit i700 wireless should not be used on patients that have been diagnosed with epilepsy due to the risk of seizures and injury. For the same reason, dental staff who have been diagnosed with epilepsy should not operate the Medit i700 wireless.

### 4.7 Explosion Hazards

#### WARNING

- The i700 wireless system is not designed to be used near flammable liquids, gases, or in environments with high oxygen concentrations.
- There is a risk of explosion if you use the i700 wireless system near flammable anesthetics.
- The rechargeable battery used with the i700 wireless is designed with safety devices included.
- The rechargeable battery shall not be exposed to excessive heat such as sunshine or the like. Failure to observe this could result in the battery explosion. Please be careful with the storage and maintenance of the battery.
- The battery charger is designed to adjust the current supply after charging is complete. But if it is not going to be used for a long time, disconnect the power to the battery charger or remove the fully charged battery from the battery charger.

### 4.8 Pacemaker and ICD Interference Risk

#### WARNING

- Do not use the i700 wireless system on patients with pacemakers and ICD devices.
- Check each manufacturer's instructions for interference by peripheral devices, such as computers used with the i700 wireless system.

### 5.1 Electromagnetic Emissions

The i700 wireless system is intended for use in the electromagnetic environment as specified below. The customer or the user of the i700 wireless system should ensure that it is used in such an environment.

| Guidance and Manufacturer's Declaration – Electromagnetic Emission |            |  |
|--|------------|--|
| Emission Test  | Compliance | Electromagnetic Environment – Guidance   |
| RF Emissions CISPR 11  | Group 1    | The i700 wireless uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.   |
| RF Emissions CISPR 11  | Class A    | The i700 wireless is suitable for use in all establishments. This includes domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic Emissions IEC 61000-3-2                                   | Class A    |  |
| Voltage Fluctuations / Flicker Emissions                           | Complies   |  |

#### WARNING

This i700 wireless system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the i700 wireless or shielding the location.

### 5.2 Electromagnetic Immunity

#### Guidance 1

The i700 wireless system is intended for use in the electromagnetic environment as specified below. The customer or the user of the i700 wireless system should ensure that it is used in such an environment.

| Guidance and Manufacturer's Declaration – Electromagnetic Immunity                                  |   |   |  |
|---|---|---|--|
| Immunity Test   | IEC 60601 Test Level  | Compliance Level  | Electromagnetic Environment – Guidance   |
| Electrostatic Discharge (ESD) IEC 61000-4-2   | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air  | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air  | Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with a synthetic material, relative humidity of at least 30% is recommended.  |
| Electrical Fast Transient / Burst IEC 61000-4-4   | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  | The mains power quality should be that of a typical commercial or hospital environment.  |
| Surge IEC 61000-4-5   | ± 0.5 kV, ± 1 kV differential mode  | ± 0.5 kV, ± 1 kV differential mode  | The mains power quality should be that of a typical commercial or hospital environment.  |
|   | ± 0.5 kV, ± 1 kV, ± 2 kV common mode  | ± 0.5 kV, ± 1 kV, ± 2 kV common mode  |  |
| Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Input Lines IEC61000-4-11 | 0% Ur (100% dip in Ur) for 0.5/1 cycles<br>70% Ur (30% dip in Ur) for 25/30 cycles<br>0% Ur (100% dip in Ur) for 250/300 cycles | 0% Ur (100% dip in Ur) for 0.5/1 cycles<br>70% Ur (30% dip in Ur) for 25/30 cycles<br>0% Ur (100% dip in Ur) for 250/300 cycles | The mains power quality should be that of a typical commercial or hospital environment. If the user of the i700 Wireless system requires continued operation during power mains interruptions, it is recommended that the i700 wireless system be powered from an uninterruptible power supply or a battery. |

|   |        |        |  |
|---|--------|--------|--|
| Power Frequency<br>Magnetic Fields<br>(50/60Hz) IEC 61000-<br>4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields<br>should be at levels characteristic of a<br>location in a typical commercial or<br>hospital environment. |
|---|--------|--------|--|

NOTE : Ur is the main voltage (AC) prior to the application of the test level.

• **Guidance 2**

| Recommended Separation Distances Between Portable and Mobile Communication Equipment and the i700 wireless |   |   |  |
|--|---|---|--|
| Rated Maximum<br>Output Power of<br>the Transmitter<br>[W]   | Separation Distance According to the Frequency of Transmitter [M] |   |  |
|  | IEC 60601-1-2: 2014   |   |  |
|  | 150 kHz to 80 MHz<br>$d = 1.2 \sqrt{P}$                           | 80 MHz to 2.7 GHz<br>$d = 2.0 \sqrt{P}$ |  |
| 0.01   | 0.12  | 0.20                                    |  |
| 0.1  | 0.38  | 0.63                                    |  |
| 1  | 1.2   | 2.0                                     |  |
| 10   | 3.8   | 6.3                                     |  |
| 100  | 12  | 20                                      |  |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 : At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

• **Guidance 3**

The i700 wireless system is intended for use in the electromagnetic environment specified below. The customer or the user of the i700 wireless system should ensure that it is used in such an environment.

| Guidance and Manufacturer's Declaration – Electromagnetic Immunity |   |                  |  |
|--|---|------------------|--|
| Immunity Test  | IEC 60601 Test Level                                  | Compliance Level | Electromagnetic Environment – Guidance   |
| Conducted RF<br>IEC 61000-4-6                                      | 3 Vrms 150 kHz to 80 MHz<br>Outside ISM Bands amateur | 3 Vrms           | Portable and mobile RF communications<br>equipment should not be used closer to any part<br>of the Ultrasound System, including cables, than<br>the recommended separation distance. This is<br>calculated using the equation applicable to the<br>frequency of the transmitter.<br><br>Recommended Separation Distance (d):<br>$d = 1.2 \sqrt{P}$<br>IEC 60601-1-2:2007<br>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.3 \sqrt{P}$ 80 MHz to 2.5 GHz<br>IEC 60601-1-2:2014<br>$d = 2.0 \sqrt{P}$ 80 MHz to 2.7 GHz                                       |
|  | 6 Vrms 150 kHz to 80 MHz<br>In ISM Bands amateur      | 6 Vrms           |  |
| Radiated RF IEC<br>61000-4-3                                       | 3 V/m 80 MHz to 2.7 GHz                               | 3 V/m            | Where $P$ is the maximum output power rating of<br>the transmitter in watts (W) according to the<br>transmitter manufacturer, $d$ is the recommended<br>separation distance in meters (m).<br>Field strengths from fixed RF transmitters, as<br>determined by an electromagnetic site survey<br>should be less than the compliance level in each<br>frequency range.<br>Interference may occur in the vicinity of<br>equipment marked with the following symbol :<br><br> |

NOTE 1 : At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

NOTE 3 : The ISM (Industrial, Scientific, and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

▪ **Guidance 4**

The i700 wireless system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the i700 wireless system. Otherwise, degradation of the performance of this equipment could result.

| Guidance and Manufacturer's Declaration – Electromagnetic Immunity |                    |   |                                   |                      |                  |
|--|--------------------|---|-----------------------------------|----------------------|------------------|
| Immunity Test  | Band <sup>1)</sup> | Service <sup>1)</sup>   | Modulation                        | IEC 60601 Test Level | Compliance Level |
| Proximity Fields from RF Wireless Communications IEC 61000-4-3     | 380 – 390 MHz      | TETRA 400   | Pulse Modulation 18 Hz            | 27 V/m               | 27 V/m           |
|  | 430 – 470 MHz      | GMRS 460; FRS 460   | FM ±5 kHz Deviation<br>1 kHz sine | 28 V/m               | 28 V/m           |
|  | 704 – 787 MHz      | LTE Band 13, 17   | Pulse Modulation 217 Hz           | 9 V/m                | 9 V/m            |
|  | 800 – 960 MHz      | GSM 800;900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5          | Pulse Modulation 18 Hz            | 28 V/m               | 28 V/m           |
|  | 1700 – 1990 MHz    | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse Modulation 217 Hz           | 28 V/m               | 28 V/m           |
|  | 2400 – 2570 MHz    | Bluetooth; WLAN 802.11b/g/n; RFID 2450; LTE Band 7              | Pulse Modulation 217 Hz           | 28 V/m               | 28 V/m           |
|  | 5100 – 5800 MHz    | WLAN 802.11a/n  | Pulse Modulation 217 Hz           | 9 V/m                | 9 V/m            |

NOTE : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

 **WARNING**

- Use of the i700 wireless adjacent to or on other equipment must be avoided as it may result in improper operation. If this use is necessary, it is advisable that this and the other equipment be observed to verify that they are operating normally.
- The use of accessories, transducers, and cables other than those specified or provided by the Medit of the i700 wireless could result in high electromagnetic emissions or reduced electromagnetic immunity from this equipment and result in improper operation.

<sup>1)</sup> For some services, only the uplink frequencies are included.

## 6. Wireless Compliance Information

### 6.1 FCC Compliance Statement

Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

**FCC No.:** 2A2QM-MD-IS0300

IMPORTANT NOTE:

FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with a minimum distance of 20 cm between the radiator and your body.

### 6.2 IC Compliance Statement

This Class B digital apparatus complies with Canadian ICES-003.

This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept

any interference, including interference that may cause undesired operation of the device. Any changes or modifications not expressly approved by the manufacturer could void the user's authority to operate this equipment.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

The device could automatically discontinue transmission in case of the absence of information to transmit or operational failure. Note that this is not intended to prohibit transmission of control or signaling information or the use of repetitive codes where required by the technology.

IC No.: 27675MD-IS0300

IMPORTANT NOTE:

IC Radiation Exposure Statement

This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 cm between the radiator and your body.

Transmit Antenna Notice

This radio transmitter (IC: 27675MD-IS0300) has been approved by Innovation, Science and Economic Development Canada to operate with the antenna types listed below, with the maximum permissible gain indicated. Antenna types not included in this list that have a gain greater than the maximum gain indicated for any type listed are strictly prohibited from use with this device.

#### Antenna List

| Model         | Type                | Maximum Gain (dBi) |
|---------------|---------------------|--------------------|
| Sil6310       | Patch Array Antenna | 18 dBi             |
| 2450AT07A0100 | Passive Antenna     | 1 dBi              |

### 6.3 KC Compliance Statement

This device has been assessed for conformity for use in business environments. When used in residential environments, there is a risk of radio wave interference.

### 6.4 TELEC (Japan) Compliance Statement

R 209-J00306  
R 209-J00282  
R 020-180117  
R 020-180116

## 7. Specifications

| Model Name  | MD-IS0300   |
|---|---|
| Trade Name  | i700 wireless   |
| Packing Unit  | 1 set   |
| Classifications for Protection Against Electric Shock | Class I, Type BF Applied Parts                          |
| * This product is a medical device.                   |   |
| Handpiece   |   |
| Dimension   | 312.7 x 43.9 x 47.4 mm (W x L x H)                      |
| Weight  | 254.5 g   |
| Rating  | 4 V <sup>~</sup> , 4 A                                  |
| DC Adapter  |   |
| Model Name  | ATM065T-P120  |
| Input   | 100 – 240 Vac, 50 – 60 Hz                               |
| Output  | 12 V <sup>~</sup> , 5 A                                 |
| Dimension   | 119 x 60 x 36 mm (W x L x H)                            |
| Rechargeable Battery                                  |   |
| Model Name  | MD-IS0300REB  |
| Type  | Lithium-ion   |
| Output  | 3.6 Vdc, 11.16 Wh                                       |
| Dimension   | 21.4 x 73.4 mm  |
| Weight  | 60 g  |
| Capacity  | 3,100 mAh   |
| Wireless Hub  |   |
| Input   | 12 V <sup>~</sup> , 5 A<br>9 V / 5 V <sup>~</sup> , 3 A |
| Dimension   | 100 x 94.8 x 44.4 mm (W x L x H)                        |
| Weight  | 181 g   |
| Battery Charger                                       |   |
| Input / Output  | 12 V <sup>~</sup> , 5 A                                 |

|   |                                  |  |
|---|----------------------------------|--|
| Dimension   | 44.7 x 100 mm (H x Ø)            |  |
| Weight  | 247 g                            |  |
| <b>Calibration Tool</b>                             |                                  |  |
| Dimension   | 123.8 x 54 mm (H x Ø)            |  |
| Weight  | 220 g                            |  |
| <b>Auto Calibration Tool (*Sold Separately)</b>     |                                  |  |
| Dimension   | 168.7 x 92.1 x 48.2 mm           |  |
| Weight  | 492 g                            |  |
| Output  | 3.6 Vdc, 11.16 Wh (MD-IS0300ECB) |  |
| <b>Wireless Module</b>                              |                                  |  |
| 60 GHz  | Frequency Bands                  | HRP: 60.48 – 62.64 GHz<br>MRP: 60.48 – 62.64 GHz<br>LRP: 60.16 – 62.96 GHz |
|   | Modulation Type                  | BPSK   |
|   | e.i.r.p.                         | HRP: 24.2 dBm<br>MRP: 24.0 dBm<br>LRP: 14.6 dBm                            |
|   | Antenna Gain                     | 18.0 dBi   |
| 2.4 GHz (Bluetooth LE)                              | Frequency Bands                  | 2402 – 2480 MHz  |
|   | Channels                         | 40   |
|   | Channel Bandwidth                | 1 MHz<br>2 MHz   |
|   | e.i.r.p.                         | 9.8 dBm<br>A-variant: 19.7 dBm<br>N-variant: 19.8 dBm                      |
|   | Modulation                       | GFSK   |
|   | Antenna Gain                     | A-variant: 1 dBi<br>N-variant: 2.14 dBi                                    |
| <b>Operating, Storage, and Transport Conditions</b> |                                  |  |
| Operating Condition                                 | Temperature                      | 18 – 28°C (64.4 – 82.4°F)  |
|   | Humidity                         | 20 – 75% relative humidity (non-condensing)                                |
|   | Air Pressure                     | 800 – 1,100 hPa  |
| Storage Condition                                   | Temperature                      | -5 – 45°C (23 – 113°F)   |
|   | Humidity                         | 20 – 80% relative humidity (non-condensing)                                |
|   | Air Pressure                     | 800 – 1,100 hPa  |
| Transport Condition                                 | Temperature                      | -5 – 45°C (23 – 113°F)   |
|   | Humidity                         | 20 – 80% relative humidity (non-condensing)                                |
|   | Air Pressure                     | 620 – 1,200 hPa  |

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