# dimao Print Denture Base

dimao Print Denture Base

- DE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- NL** Gebruiksaanwijzing

Explanation of symbols on labelling	
	Batch code
	Manufacturer
	Keep away from sunlight
	Consult instructions for use
	Use-by date
	Storage temperature

Manufacturer:  
DENTCA, INC.  
357 VAN NESS WAY, STE 250  
TORRANCE, CA 90501  
USA

Distributor:  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)  
+49 800 4372522

Distributed in USA /  
Canada exclusively by:  
Kulzer, LLC  
4315 South Lafayette Blvd.  
South Bend, IN 46614-2517  
1-800-431-1785

Caution: Federal law restricts  
this device to sale by or on the  
order of a dental professional.

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EU Representative:  
MT Promed Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert  
Germany

Medizinprodukt  
Medical device  
Dispositivo médico  
Producto sanitario  
Dispositivo medico  
Equipamento para saúde  
Medisch hulpmiddel

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## GB Instructions for use dimao Print Denture Base

dimao Print Denture Base is a light-curable resin indicated for fabrication and repair of full and partial removable dentures and basesets. The material is an alternative to traditional heat-curable and auto polymerizing resins. Fabrication of dental prosthetics with dimao Print Denture Base requires a computer-aided design and manufacturing (CAD/CAM) system that includes the following components: digital denture base files based on a digital impression, stereolithographic additive printer, and curing light equipment.

- Requirements**
1. Digital denture base file; STL format file
  2. Additive printer and its operation software;

	Additive Printer	Operation Software	Provider
	cara Print	cara Print CAM	Kulzer

	Curing Equipment	Provider
	HiLite power 3D	Kulzer

### Specific Manufacturing Considerations

1. Digital denture base file
  - 1.1 File format: STL file
  - 1.2 Digital design: denture base or baseslates with following minimum thickness of the areas:

Area	Lingual Ridge area	Palatal / Lingual	Fazial / Buccal	IOD* area
Maxillary	≥ 2,5 mm	≥ 2 mm (≥ 3 mm single Arch upper)	≥ 2 mm	≥ 2,5 mm
Mandibular	≥ 2,5 mm	≥ 2 mm	≥ 2 mm	≥ 2,5 mm

- 1.3 File size: upload-able in the 3D printer operation software

2. Additive printer
  - 2.1 Hardware
    - a. Wellenlänge: 385 nm or 405 nm
    - b. Light source
      - Stereolithographische (SLA) method; laser with 25 mW < X < 250 mW
      - Digital Light Processing (DLP) method
    - c. Build Volume: > 103 x 58 x 130 mm (Least fit one arch)
    - d. Build Path: line drawing path or surface drawing path
  - 2.2 Features of Operation Software
    - a. STL file import
    - b. Automatic rotation and placement
    - c. Layer slicer for path inspection
    - d. Auto and manual generation of supports
  - 2.3 Printing Parameters

Printer Model	Layer Thickness (micron)	Recommended orientation angle (degree)	Support point size (mm)	Support density (mm)
cara Print 4.0	50, 100	20–90	2–4	min 1.5

- 2.4 Environmental Conditions
  - a. Temperature: 18–30°C (64–86°F)
  - b. Relative Humidity: 30–90%
- 2.5 Cleaning Kit
  - Rinse bath and tubs, paper towel, squeeze bottle for isopropyl alcohol, Scraper
3. Recommended Curing Light Equipment (Post curing units)
  - 3.1 Flood Type Curing Equipment

Manufacturer / Model	Curing Chamber	Supply voltage (V / Hz)	Lamp power (W)	Light intensity (mW/cm²)	Lamp wavelength (nm)	Curing time (min)
Kulzer / HiLite power 3D	Required	100,115,230 / 50–60	200	N/A	390–540	20

- 3.2 Accessories
  - a. USP Grade glycerin
  - b. Transparent glass container and 2 glass plates
  - c. Heat-protective gloves and silicone coated stainless steel tong
  - d. Thermocouple
4. Notification  
The device specifications have been validated using the software, printers, and process parameters specified in this document. Any other printers, operation software and post-printing processes will be outside of the device specifications and the FDA clearance. Users shall follow this document to use the device.

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- Warnings:**
1. dimao Print Denture Base contains polymerizable monomers which may cause skin irritation (allergic contact dermatitis) or other allergic reactions in susceptible persons. If contact with skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
  2. Avoid inhalation or ingestion. High vapor concentration can cause headache, irritation of eyes or respiratory system. Direct contact with eyes may cause possible corneal damage. Long-term excessive exposure to the material may cause serious health effects. Monitor air quality per OSHA standards. Eye Contact: Immediately flush eyes with plenty of clean water for at least 20 minutes, and consult a physician. Wash the contacted area thoroughly with soap and water; irrigation: in case of exposure to a high concentration of vapor or mist, remove person to fresh air. Give oxygen or artificial respiration as required. Ingestion: Contact your regional poison control center immediately.
  3. BURN HAZARD: Glycerol Bath CAN REACH TEMPERATURES OF 90°C (~200°F) AND LEAD TO SEVERE BURNS. Only trained users should perform the glycerol curing step with caution and appropriate PPE. We also recommend placing a warning label on the window of the cure unit to alert all lab users to the potential hazard.

- Precautions:**
1. When washing the printed denture base with solvent or grinding the denture base, it should be in a properly ventilated environment with proper protective masks and gloves.
  2. Store dimao Print Denture Base at or below 15–25°C (60–77°F) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date.
  3. Exposed or unused dimao Print Denture Base should be completely cured or polymerized prior disposal.

- Adverse Reactions:**
1. Direct contact with the uncured resin may induce skin sensitization in susceptible individuals.
  2. Proper ventilation and personal protective equipment should be used when grinding denture base resins as the particulate generated during grinding may cause respiratory, skin and eye irritation.

### Procedure to Fabricate Denture

1. Printing Preparation
  - a. Select the denture base shade based on prescription. (Recommended to use the different resin tank for the different shade.)
  - b. Open the 3D printer cover and fill the resin tank of the printer with fresh dimao Print Denture Base up to the required filling line by manufacturer. (When filling the resin into the resin tank, gloves and mask should be used.)
  - c. Close the printer cover.
2. Printing
  - a. Load the denture base model file to be printed in printer operation software which printer manufacturer recommended.
  - b. Use auto-orientation or manual orientation to find its optimal position for printing. The recommended orientation by printer provider is a tilted orientation such as space diagonal from 20 to 90 angle. If auto-orientation is not satisfied, rotate to make optimal position.
  - c. Generate support structures on the denture base using the recommended setting by printer provider such as support contact size and support density. (Avoid the support structures inside the tooth socket area and the ridge area. If some of supports are generated inside the tooth sockets and on the ridge area, either relocate it to the edge of denture base or remove it.)
  - d. Use layout tools of the software to move the denture base model within the imaginary build platform to prevent the overlapping between models.
  - e. Start printing.
3. Cleaning
  - a. Detach the printed denture base from the build platform.
  - b. Use a small flush cutter to remove the support sticks from the denture base.
  - c. Wash the denture base with isopropyl alcohol.
  - d. Use a paper towel to dry the denture base or dry it at room temperature under ventilation system or open area.

### Denture Fabrication using a printed Denture Teeth and base

1. Bonding the printed Teeth to the printed denture base
  - a. Prepare, before post-cured, the printed Teeth and printed denture base with socket shapes to receive printed Teeth (Tooth Sockets).
  - b. Place the printed Teeth onto the dimao Print Denture Teeth into the corresponding Tooth Sockets on the printed denture base and check Teeth fitting.
  - c. Apply the small amount of light curable adhesive into the Tooth Sockets and bond Teeth by exposing into UV light until the Teeth set in position. Recommended to use dimao Print Denture Base for bonding.
  - d. If necessary, apply small amount of dimao Print Denture Base using an applicator to smoothen the edges of the denture base and cure it.
2. Post Curing
  - a. Cure the final denture by sinking into the glycerin container (Glycerin temperature should be greater than 60°C (140°F) and it is recommended the glycerine to replace every 80 hrs running or every three months whichever comes first), for the required curing time under recommended post-curing unit. For all of the post-curing time tissue sub to and for another half of the time tissue sub down.
  - b. Take out the printed denture from the curing oven using coated tong (Be careful hot glycerin).
  - c. Rinse the cured denture with a water.
3. Finishing
  - a. Smoothen the support spots on the denture base using a bur.
  - b. Polish the final denture with wet polishing sand by conventional method.

### Repair of denture and baseslates

- Repair:** This process is only applied to temporary repair cases. It is recommended to remake the whole denture using an original design file.
1. Prepare a cast made of a putty using a broken denture.
  2. Prepare the fracture area by grinding to open more and roughing the outer side of fracture area.
  3. Prime the roughened surfaces of the repair area with dimao Print Denture Base.
  4. Place the broken denture on the cast.
  5. Apply the dimao Print Denture Base resin to cover the roughed and fracture areas and cure the areas by exposing light curing machine until the resin is solidified.
  6. Place the denture on the cast in the post-curing machine for the half of required time. Carefully remove the denture from the cast and cure the tissue side for the half of required time.
  7. Grind, polish and finish.

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