

APPLICATION GUIDE

FOR DIMA PRINT STONE GRAY

Version 06.2022



MORE INFO
KULZER.COM/
STONEGRAY



KULZER
MITSUI CHEMICALS GROUP

Equipment you need

3shape
exocad

Design software

CAD design software (3shape) and exocad



Print resin for 3D printing

dima Print Stone gray



3D printer

cara Print 4.0 pro



Wash unit

cara Print Clean pro or ultra sonic bath



Post curing unit

cara Print LEDcure or HiLite power 3D

1. Design notes

Recommended settings for exocad



Plateless model design					
Dies	Pin height	Preparation margin	Seating width	Shaft taper	Taper height limit
	2.5mm	0mm	1.5mm	0°	10mm
	<input checked="" type="checkbox"/> Add concavity		<input checked="" type="checkbox"/> Die bottom parallel to model base		
Settings	Base				
	Horizontal shaft gap		Vertical shaft gap		Pedestal height
	-0.015mm		-0.015mm		0mm
	<input checked="" type="checkbox"/> Hollow model				
	Wall thickness		Bottom sill		Cavity fill diameter
3mm		1mm		0mm	


Recommended settings for 3Shape



Model recommended wall thickness	Partial model minimum wall thickness
3-4mm	2mm

2. CAM design for cara Print 4.0 pro with cara CAM 2.0

Print position [angle and orientation]



The recommended orientation is flat on the build table with 0° angle for all model types.

2.a Die models

SUPPORTS: no supports needed.

Supports




Top view



Bottom view


! NOTE:
NO SUPPORTS NEEDED IN THE YELLOW AREAS.



! NOTE:
THE BEST PRINTING RESULT FOR 2 MODELS CAN BE ACHIEVED BY POSITIONING THE U-SHAPED ARCHES AS IN THE PICTURE.


2.b Full model – hollow

SUPPORTS: it is important to manually add supports to the inside of hollow models.



Right click to add supports

SUPPORTS: Set the supports automatically via a right click (in the yellow marked area).



Section view: with supports in the model

3. Processing

Mixing of photopolymer

- ▶ Time: 5 min
- ▶ Shake well before pouring the liquid into printer vat
- ▶ Liquid is quite viscous, all contained pigments needed to be shaken up otherwise printing errors can occur

Printing parameters

- ▶ Choose printing parameters for dima Print Stone gray
- ▶ Resolution: 50 µm/70 µm/100 µm

Printing performance overview


Example	Build height [mm]	Time [min] for 50 µm resolution	Time [min] for 70 µm resolution	Time [min] for 100 µm resolution	Material consumption [ml]
One upper jaw full (hollow)	28.7	66	44	30	29
One lower jaw full (hollow)	27.9	61	41	29	22.5
One upper jaw with 6 dies	23.2	55	37	26	27
One lower jaw with 7 dies	24.9	62	41	28	20.8

4. Cleaning and post-curing

Cleaning recommendations

Please use 99% IPA to clean dima Print Stone gray in a cleaning unit or an ultrasonic bath.

We generally recommend using cara Print Clean pro to wash printed objects. For more information, see the instructions for the cleaning unit on our website: kulzer.com/cleanpro.



Step	Cleaning	Curing
1.	 <p>cara Print Clean pro (or ultra sonic bath)</p>	 <p>OR</p>  <p>HiLite Power 3D</p>
	<p>1. Pre-cleaning: 3 min. 2. Post-cleaning: 2 min.</p>	<p>Select dima Print Stone gray program (no turning of object necessary)</p> <p>10 minutes (5 minutes front side + 5 minutes reverse side)</p>

PERFECT MATCH: DIMA PRINT STONE GRAY & HERACERAM SAPHIR

dima Print materials are light-curing monomers specially optimized for reliable results using the 3D-Printing pro Solution.

TESTED BY DENTAL
EXPERTS



BROAD VARIETY
OF INDICATIONS

IDEAL PROPERTIES
FOR EACH INDICATION
VISIT
[WWW.KULZER.COM/
DIMA](http://WWW.KULZER.COM/DIMA)
FOR A COMPLETE
OVERVIEW.



KULZER
MITSUI CHEMICALS GROUP



LEARN MORE
ABOUT
DIMA PRINT
STONE GRAY
[KULZER.COM/
STONEGRAY](http://KULZER.COM/STONEGRAY)

Contact in Germany
Kulzer GmbH
Leipziger Straße 2
63450 Hanau, Germany