dima Mill Zirconia Disc

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dima Mill Zirconia discs

KULZER
MITSUBISHI CHEMICALS GROUP

Manufacturer:
Metoxit AG
Emdwiesenstrasse 6
8240 Thayngen (Switzerland)

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Instructions for use

Indications

dima Zirconia discs are indicated for use as a substructure for porcelain fused ceramic fixed dental restorations. Limitations are listed in table 1.

Table 1: Indications of use and maximum number of pontics.

Process chain	Single unit crowns		Bridges (maximum number of pontics in between)	
	Anterior	Posterior	Anterior	Posterior
dima Mill Zirconia ST	X	X	2	2
dima Mill Zirconia HT/HT+/ML+	X	X	2	2
dima Mill Zirconia HTE*	X	X	1	1

* dima Mill Zirconia HTE is released for crowns and bridges with max. 3 units.

Contraindications

Contraindications are all applications not indicated in table 1 above.

Application

1. Fix the disc with the printed side (occlusal) upwards according to the CAM-System instructions, set the enlargement-factor.
2. Start the machining.

3. After machining, dismount the blank with the framework.
4. Remove the framework from the blank with a diamond milling cutter.
5. Remove dust the framework gently, with a soft brush or an air-gun.***
6. Position the dry framework in a sinter-bed inside the crucible.
7. Programme your furnace according to its manual and dima sinter instructions.
8. Sinter the framework.
9. Wait until cool down is complete before removing the sintered framework.
10. Inspect the framework in respect of flaws, wall thickness, seat and margin fit. If necessary, make small adjustments e.g. with a water-cooled diamond-tool or another Zirconium dioxide suitable tool.
11. Clean the framework with water and dry it. It is now ready for veneering!

Sinter instructions

dima Mill Zirconia (ST/HT/HT+/ML+/HTE)	Ramp Rate °C/min	Final temperature T _{max} °C	Dwell time min	Cooling °C/min
Standard	5	1450–1530	120	10
Economic*	10	1450–1530	120	10 (20)***
Speed**	45	1530	30	45

* Bridges with max. 7 units

** Bridges with max. 3 units

*** In sintering furnaces with two-stage cooling programming, cool down at a rate of 10 °C/min above 900 °C, and cool down at a rate of 20 °C/min below 900 °C.

Please note as well the sintering guideline on www.kulzer.com/dima in the download area.

Safety instructions

- *** Processing dima Mill Zirconia Disc generates dust that can irritate the eyes, skin and airways. Therefore always ensure that the suction is working correctly.
- Design the framework to maximise wall thickness and cross-sections.
- Carefully read the MSDS.

Dated: 2021-06