## dima Mill CoCr solid - Rohlinge

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## dima CoCr solid blanks

**GB** Instructions for use

**Indications:** dima Mill CoCr solid blanks are indicated as Type 4 for the manufacture of fixed and removable dentures for crowns and bridge frameworks as well as bridge and implant superstructures according to DIN EN ISO 22674 and EN ISO 9693. The processing is carried out using CAD/CAM technology with suitable milling machines. (Please note manufacturer's specifications).

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X = approved

**Contraindications:** All applications not listed under Indications or in Table 1 are contraindicated.

In the case of known hypersensitivity (allergy) to the constituents of the alloy.

**Application**

**1. Fabrication:** Insert the disc into the machine according to the operating instructions of the CAM system and process the disc using tungsten carbide burs suitable for CoCr. When finished, remove the blank with the frameworks produced.

**2. Separating the frameworks from the blank:** The milled restorations can be cut out of the blank using cross-cut tungsten carbide burs or separating discs suitable for CoCr alloys.

**3. Finishing / cleaning:** Dress the separated frameworks with tungsten carbide burs being careful not to apply excess pressure. Doing so may cause surface flow of the alloy. Trim in one direction only and avoid the creation of sharp edges. For the desired veneering, the surfaces must be blasted with aluminum oxide (110-125 µm) at a pressure of 3-4 bar, thoroughly cleaned and then degreased with ethanol. △ Metal dust is harmful to health. During finishing procedures, use suitable suction and a respirator mask of type FFP3-EN 149-2001.

**4. Oxidizing:** An oxide firing is not necessary with dima Mill CoCr solid. If, however, an oxide firing (rapid heating to 1000 °C without vacuum, cooling down immediately after reaching the temperature) is carried out for visual inspection of the framework conditioning, the oxide must subsequently be stripped off with single-use abrasive aluminum oxide (110-125 µm) and subsequently cleaned (e.g. by evaporation).

**5. Joining (soldering/lasering):** For soldering, a solder and flux suitable for CoCr alloys such as Heraenium P Lot 1 in combination with Hera SLP 99 flux is recommended. Keep the soldering block as small as possible and heat it initially at about 600 °C for 10 minutes in the oven. Soldering gaps should be between 0,05 and 0,2 mm to avoid solder flowing through solder joints. For laser welding prior to firing, a special laser welding wire such as Kulzer CoCr laser welding wire (diameter 0,5 mm or 0,35 mm) is recommended.

**6. Veneering with ceramics:** We recommend the use of e.g. HeraCeram. Follow all procedures as indicated in the instruction booklet. If other ceramics are used, the instructions for use from those ceramic manufacturers must be followed.

**7. Polishing:** After ceramic firing or soldering, existing oxide/flux residues must be removed and the framework finished to a high gloss using suitable grinding/polishing instruments and polishing pastes.

**Technical Specifications:**

Alloy type	4
Density	8,8 g/cm³
CTE value (25-500 °C)	14,4*10⁻⁶ K⁻¹
Mass content of all components (in wt%)	Co 63, Cr 29, Mo 6, (Mn, Nb, Si, Fe) <1
Vickers hardness	350 HV10
Tensile strength	680 MPa
0,2 % yield strength	400 MPa
Elongation	9%
Young's modulus of elasticity	250 GPa

**Side effects:** In individual cases, hypersensitivity reactions to constituents of the alloy used cannot be ruled out. The product should not be used in case of known allergies to the constituents of the alloy.

**Interactions:** On contact with or in close proximity to alloys of a different type, galvanic effects can sometimes occur. In the event of sustained local discomfort, other materials must be used for the denture.

**Safety instructions**

• When designing the prosthetic dentures, ensure adequate wall thickness and connector cross-sections.

• These medical devices may only be used and applied by qualified staff in a dental laboratory.

Please also observe the material safety data sheet.

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