User Manual

i600

Revision 1 (January 2022)

CONTENTS

| | |
|-----|-----------|
| 02 | English |
| 24 | العربية |
| 46 | 中文(简体) |
| 66 | 中文(繁體) |
| 86 | Dutch |
| 110 | Français |
| 134 | Deutsch |
| 158 | Italiano |
| 182 | 日本語 |
| 204 | Português |
| 228 | Polskie |
| 250 | русский |
| 274 | Español |
| 298 | Svenska |
| 320 | 한국어 |

English

| | | | | | |
|----------|---|-----------|----------|---|-----------|
| 1 | About this guide | 4 | 5 | Safety Guide | 14 |
| 2 | Introduction and Overview | 4 | 5.1 | System Basics | 14 |
| 2.1 | Intended Use | 4 | 5.2 | Proper Training | 15 |
| 2.2 | Indication for Use | 4 | 5.3 | In Case of Equipment Failure | 16 |
| 2.3 | Contraindications | 5 | 5.4 | Hygiene | 16 |
| 2.4 | Qualifications of the Operating User | 5 | 5.5 | Electrical Safety | 16 |
| 2.5 | Symbols | 5 | 5.6 | Eye Safety | 17 |
| 2.6 | i600 Components Overview | 6 | 5.7 | Explosion Hazards | 17 |
| 2.7 | Setting up the i600 Device | 7 | 5.8 | Pacemaker and ICD Interference Risk | 17 |
| 2.7.1 | Basic settings of i600 | 7 | 6 | Electro – Magnetic Compatibility Information | 18 |
| 2.7.2 | Placing on Desktop Cradle | 8 | 6.1 | Electro-Magnetic Emissions | 18 |
| 2.7.3 | Installation of Wall Mount Holder | 9 | 6.2 | Electro-Magnetic Immunity | 18 |
| 3 | Image Acquisition Software Overview | 9 | 7 | Specification | 22 |
| 3.1 | Introduction | 9 | | | |
| 3.2 | Installation | 9 | | | |
| 3.2.1 | System Requirements | 9 | | | |
| 3.2.2 | Installation Guide | 10 | | | |
| 4 | Maintenance | 11 | | | |
| 4.1 | Calibration | 11 | | | |
| 4.2 | Cleaning, Disinfection, Sterilization Procedure | 12 | | | |
| 4.2.1 | Reusable tip | 12 | | | |
| 4.2.2 | Disinfection and Sterilization | 12 | | | |
| 4.2.3 | Mirror | 13 | | | |
| 4.2.4 | Handpiece | 13 | | | |
| 4.2.5 | Other Components | 13 | | | |
| 4.3 | Disposal | 14 | | | |
| 4.4 | Updates in Image Acquisition Software | 14 | | | |

1 About this guide

Convention in this guide

This user guide uses various symbols to highlight important information so as to ensure correct usage, prevent injury to the user and others, and prevent property damage. The meanings of the symbols used are described below.

WARNING

The WARNING symbol indicates information that, if ignored, could result in a medium risk of personal injury.

CAUTION

The CAUTION symbol indicates safety information that, if ignored, could result in a slight risk of personal injury, property damage, or damage to the system.

TIPS

The TIPS symbol indicates hints, tips, and additional information for optimal operation of the system.

2 Introduction and Overview

2.1 Intended Use

The i600 system is a dental 3D scanner intended to be used to digitally record topographical characteristics of teeth and surrounding tissues. The i600 system produces 3D scans for use in computer-assisted design and manufacturing of dental restorations.

2.2 Indication for Use

The i600 system should be used on patients who require 3D scanning for dental treatments such as:

- Single custom abutment
- Inlays & Onlays
- Single Crown
- Veneer
- 3 Unit Implant Bridge
- Up to 5 Unit Bridge
- Orthodontics
- Implant Guide
- Diagnosis Model

The i600 system can also be used in full-arch scans but various factors (intraoral environment, operator's expertise, and laboratory workflow) may affect the final results.


2.3 Contraindications


- The i600 system is not intended to be used to create images of the internal structure of teeth or the supporting skeletal structure.
- It is not intended to be used for cases with more than (4) subsequent edentulous tooth positions.

2.4 Qualifications of the Operating User

- The i600 system is designed for use by individuals with professional knowledge in dentistry and dental laboratory technology.
- The user of the i600 system is solely responsible for determining whether or not this device is suitable for a particular patient's case and circumstances.
- The user is solely responsible for the accuracy, completeness, and adequacy of all data entered into the i600 system and the provided software. The user should check the accuracy of the results and assess each individual case.
- The i600 system must be used in accordance with its accompanying user guide.
- Improper use or handling of the i600 system will void its warranty. If you require additional information on the proper use of the i600 system, please contact your local distributor.
- The user is not allowed to modify the i600 system.

2.5 Symbols

| No | Symbol | Description |
|----|---|---|
| 1 |  | The serial number of the object |
| 2 |  | Date of manufacture |
| 3 |  | Manufacturer |
| 4 |  | Caution |
| 5 |  | Warning |
| 6 |  | Instructions for User Manual |
| 7 |  | The official mark of Europe Certificate |
| 8 |  | Authorized Representative in the European Community |
| 9 |  | BF Type of applied part |
| 10 |  | WEEE Mark |
| 11 |  | Prescription use (U.S.A) |
| 12 |  | MET mark |
| 13 |  | AC |

| | | |
|----|---|---------------------------------|
| 14 |  | DC |
| 15 |  | Protective Earth (ground) |
| 16 |  | Temperature Limitation |
| 17 |  | Humidity Limitation |
| 18 |  | Atmospheric Pressure Limitation |
| 19 |  | Fragile |
| 20 |  | Keep Dry |
| 21 |  | Position |
| 22 |  | Seven-layer stacking prohibited |

2.6 i600 Components Overview

| No | Item | Qty | Appearance |
|----|----------------|-----|---|
| 1 | i600 Handpiece | 1ea |  |
| 2 | Power Hub | 1ea |  |

| | | | |
|----|----------------------|-----|---|
| 3 | i600 Handpiece Cover | 1ea |  |
| 4 | Reusable Tip | 4ea |  |
| 5 | Calibration Tool | 1ea |  |
| 6 | Practice Model | 1ea |  |
| 7 | Wrist Strap | 1ea |  |
| 8 | Desktop Cradle | 1ea |  |
| 9 | Wall Mount Holder | 1ea |  |
| 10 | Power Delivery Cable | 1ea |  |

| | | | |
|----|---|-----|---|
| 11 | USB 3.0 Cable | 1ea |  |
| 12 | Medical Adapter | 1ea |  |
| 13 | Power Cord | 1ea |  |
| 14 | USB Flash Drive (Pre-loaded with image acquisition software) | 1ea |  |
| 15 | User Guide | 1ea |  |

2.7 Setting up the i600 Device

2.7.1 Basic settings of i600

① Connect the USB C Cable to the Power Hub.

② Connect Medical Adapter to the Power Hub.

③ Connect Power Cord to the Medical Adapter.

④ Connect Power Cord to a Power Source.

⑤ Connect the USB C Cable to a PC.

Turn on the i600

Press the power button on the i600.

Wait until the USB connection indicator turns blue.

Turn on the i600

Press and hold the power button of the i600 for 3 seconds

2.7.2 Placing on Desktop Cradle.

2.7.3 Installation of Wall Mount Holder.

3 Image Acquisition Software Overview

3.1 Introduction

The image acquisition software provides a user-friendly working interface to digitally record topographical characteristics of teeth and surrounding tissues using the i600 system.

3.2 Installation

3.2.1 System Requirements

Minimum System Requirements

| | Laptop | Desktop |
|---------|--|--|
| CPU | Intel Core i7 - 10750H | Intel Core i7 - 10700K |
| | AMD Ryzen 7 4800H | AMD Ryzen 7 3800X |
| RAM | 32 GB | 32 GB |
| Graphic | Nvidia GeForce RTX 1660/2060/3060 Above 6GB (Not supporting Radeon) | Nvidia GeForce RTX 1660/2060/3060 Above 6GB (Not supporting Radeon) |
| | | |
| OS | Windows 10 Pro 64-bit / Windows 11 Pro 64-bit | |

Recommend System Requirements

| | Laptop | Desktop |
|---------|--|--|
| CPU | Intel Core i7 - 11800H | Intel Core i7 - 11700K |
| | AMD Ryzen 7 5800H / AMD Ryzen 9 5900H | AMD Ryzen 7 5800X |
| RAM | 32 GB | 32 GB |
| Graphic | Nvidia GeForce RTX 2070/2080/3070/3080/3090 Above 8GB (Not supporting Radeon) | Nvidia GeForce RTX 2070/2080/3070/3080/3090 Above 8GB (Not supporting Radeon) |
| | | |
| OS | Windows 10 Pro 64-bit / Windows 11 Pro 64-bit | |

 Use PC and monitor certified IEC 60950, IEC 55032, IEC 55024


 The USB 3.0 Cable provided with the i600 is a special cable that provides Power Delivery. With the Power Delivery System, power

can be supplied without using the supplied Power Hub, so you can scan.

When using cables other than the USB 3.0 Cable provided by MEDIT the device may not work. MEDIT is not responsible for any problems caused by cables other than the USB 3.0 Cable provided by MEDIT. Be sure to use only the USB 3.0 cable included in the package.

3.2.2 Installation Guide

- ① Run Medit_Scan_for_Clinics_x.x.x.exe

- ② Select the setup language and click "Next."

- ③ Select the installation path.

- ⑤ It may take several minutes to finish the installation process. Please do not shut down the PC until the installation is complete.

- ④ Read the "License Agreement" carefully before checking "I agree to the License terms and conditions," and then click Install.

- ⑥ After the installation is complete, restart the PC to ensure optimal program operation.

If the scanner is connected, please disconnect the scanner from the PC by removing the USB cable.

4 Maintenance

CAUTION

- Equipment maintenance should only be carried out by a MEDIT employee or a MEDIT-certified company or personnel.
- In general, users are not required to perform maintenance work on the i600 system besides calibration, cleaning, and sterilization. Preventive inspections and other regular maintenance are not required.

4.1 Calibration

Periodic calibration is required to produce precise 3D models.

You should perform a calibration when:

- The quality of the 3D model is not reliable or accurate when compared to previous results.
- Environmental conditions such as temperature have changed.

- Calibration period has expired. You can set the calibration period in the Menu > Settings > Calibration Period (Days).
- The calibration panel is a delicate component. Do not directly touch the panel. Check the calibration panel if the calibration process is not performed properly. If the calibration panel is contaminated, please contact your service provider.
- We recommend performing a calibration periodically. You can set the calibration period via Menu > Settings > Calibration Period (Days). The default calibration period is 14 days.

How to calibrate i600

- Turn on the i600 and launch the image acquisition software.
- Run the Calibration Wizard from Menu > Settings> Calibration.
- Prepare the Calibration Tool and the i600 handpiece.
- Turn the dial of the calibration tool to position **1**.
- Put the handpiece into the calibration tool.
- Click “Next” to start the calibration process.
- When the calibration tool is mounted properly in the correct position, the system will automatically acquire the data. **1**.
- When data acquisition is complete at position **1**, turn the dial to the next position.
- Repeat the steps for positions **2** ~ **8** and the **LAST** position.
- When data acquisition is complete at the **LAST** position, the system will automatically calculate and show the calibration results.

4.2 Cleaning, Disinfection, Sterilization Procedure

4.2.1 Reusable tip

The reusable tip is the part which is inserted into the patient's mouth during scanning. The tip is reusable for a limited number of times. The tip needs to be cleaned and sterilized between patients to avoid cross-contamination.

- The tip should be cleaned manually using a disinfecting solution. After cleaning and disinfecting, inspect the mirror inside the tip to ensure that there are no stains or smudges.
- Repeat the cleaning and disinfection process if necessary. Carefully dry the mirror using a paper towel.
- Insert the tip into a paper sterilization pouch and seal it, making sure that it is airtight. Use either a self-adhesive or heat-sealed pouch.
- Sterilize the wrapped tip in an autoclave with the following conditions:
 - » Sterilize for 30 minutes at 121°C (249.8°F) at gravity type and dry for 15 minutes.
 - » Sterilize for 4 minutes at 134°C (273.2°F) at pre-vacuum type and dry for 20 minutes.
- Use an autoclave program that dries the wrapped tip before opening the autoclave.
- Scanner tips can be re-sterilized up to 100 times and must thereafter be disposed of as described in the disposal section.

4.2.2 Disinfection and Sterilization

- Clean the tip immediately after use with soap water and a brush. We recommend using a mild dishwashing liquid. Make sure the mirror of the tip is completely clean and stain-free after cleaning. If the mirror appears stained or foggy, repeat the cleaning process and rinse thoroughly with water. Dry the mirror carefully with a paper towel.
- Disinfect the tip using Wavicide-01 for 45 to 60 minutes. Please refer to the Wavicide-01 Solution instruction manual for proper use.
- After 45 to 60 minutes, remove the tip from the disinfectant and rinse thoroughly.
- Use a sterilized and non-abrasive cloth to gently dry the mirror and the tip.

CAUTION

- The mirror found in the tip is a delicate optical component which should be handled with care to ensure optimal scan quality. Be careful not to scratch or smudge it as any damage or blemishes may affect the data acquired.
- Make sure to always wrap the tip before autoclaving. If you autoclave an exposed tip, it will cause stains on the mirror which cannot be removed. Check the autoclave manual for more information.
- New tips need to be cleaned and sterilized / autoclaved before their first use.
- Medit is not responsible for any damage including distortion,

blackening etc.

4.2.3 Mirror

The presence of impurities or smudges on the tip mirror may lead to poor scan quality and an overall poor scanning experience. In such a situation, you should clean the mirror following the steps below:

- Disconnect the scanner tip from the i600 handpiece.
- Pour alcohol on a clean cloth or cotton-tipped swab and wipe the mirror. Make sure to use alcohol that is free of impurities or it may stain the mirror. You can use either ethanol or propanol (ethyl-/propyl alcohol).
- Wipe the mirror dry using a dry, lint-free cloth.
- Make sure the mirror is free of dust and fibers. Repeat the cleaning process as necessary.

4.2.4 Handpiece

After treatment, clean and disinfect all other surfaces of the handpiece except for the scanner front (optical window) and end (air vent hole).

Cleaning and disinfecting must be done with the device turned off. Use the device only after it is completely dry.

Recommended cleaning and disinfecting solution:

Denatured alcohol (aka. ethyl alcohol or ethanol) – typically 60-70% Alc/Vol.

The general cleaning and disinfecting procedures are as follows:

- Turn off the device using the power button.
-

-
- Unplug all the cables from the power hub.
 - Attach the handpiece cover to the front of the scanner.
 - Pour the disinfectant onto a soft, lint-free, and non-abrasive cloth.
 - Wipe the scanner surface with the cloth.
 - Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

CAUTION

- Do not clean the handpiece when the device is turned on as the fluid may enter the scanner and cause malfunction.
- Use the device after it is completely dry.

CAUTION

- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

4.2.5 Other Components

- Pour the cleaning and disinfecting solution onto a soft, lint-free, and non-abrasive cloth.
- Wipe the component surface with the cloth.
- Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

CAUTION

- Chemical cracks may appear if an improper cleaning solution is used during cleaning.
-

4.3 Disposal

CAUTION

- The scanner tip must be sterilized before disposal. Sterilize the tip as described in section 4.2.1.
- Dispose of the scanner tip as you would any other clinical waste.
- Other components are designed to conform with the following directives:
- RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment. (2011/65/EU)
- WEEE, Waste Electrical and Electronic Equipment Directive. (2012/19/EU)

4.4 Updates in Image Acquisition Software

The image acquisition software automatically checks for updates when the software is in operation.

If a new version of the software is released, the system will automatically download it.

5 Safety Guide

Please adhere to all the safety procedures as detailed in this user guide to prevent human injury and equipment damage. This document uses the words WARNING and CAUTION when highlighting precautionary messages.

Carefully read and understand the guidelines, including all precautionary messages as prefaced by the words WARNING and CAUTION. To avoid bodily injury or equipment damage, make sure to adhere strictly to the safety guidelines. All instructions and precautions as specified in the Safety Guide must be observed to ensure proper functionality of the system and personal safety.

The i600 system should only be operated by dental professionals and technicians who are trained to use the system. Using the i600 system for any purpose other than its intended usage as outlined in section "2.1 Intended Use" may result in injury or damage to the equipment. Please handle the i600 system according to the guidelines in the safety guide.

5.1 System Basics

CAUTION

- The USB 3.0 cable connected to the Power Hub is the same as a standard USB cable connector. However, the device may not operate normally if a standard 3.0 USB cable is used with i600.
 - The connector provided with the Power Hub is designed specifically for the i600 and should not be used with any other device.
 - If the product has been stored in a cold environment, give it time to adjust to the temperature of the environment before use. If used immediately, condensation may occur which may damage the electronic parts inside the unit.
 - Ensure that all components provided are free from physical damage. Safety cannot be guaranteed if there is any physical
-

damage to the unit.

- Before using the system, check that there are no issues such as physical damage or loose parts. If there is any visible damage, do not use the product and contact the manufacturer or your local representative.
 - Check the i600 body and its accessories for any sharp edges.
 - When not in use, the i600 should be kept mounted on a desk stand or wall mount stand.
 - Do not install the desk stand on an inclined surface.
 - Do not place any object on the i600 body.
 - Do not place the i600 on any heated or wet surface.
 - Do not block the air vents located at the rear of the i600 system. If the equipment overheats, the i600 system may malfunction or stop working.
 - Do not spill any liquid on the i600 device.
 - Do not pull or bend the cable connected to the i600.
 - Carefully arrange all the cables so that you or your patient do not trip or get caught in the cables. Any pulling tension on the cables may cause damage to the i600 system.
 - Always place the power cord of the i600 system in an easily accessible location.
 - Always keep an eye on the product and your patient while using the product to check for abnormalities.
 - If you drop the i600 tip on the floor, do not attempt to reuse it. Discard the tip immediately as there is a risk that the mirror attached to the tip may have been dislodged.
-

- Due to its fragile nature, the i600 tips should be handled with care. To prevent damage to the tip and its internal mirror, be careful to avoid contact with a patient's teeth or restorations.
- If the i600 is dropped on the floor or if the unit is impacted, it must be calibrated before use. If the instrument is unable to connect to the software, consult the manufacturer or authorized resellers.
- If the equipment fails to operate normally, such as having issues with accuracy, stop using the product, and contact the manufacturer or authorized resellers.
- Install and use only approved programs to ensure proper functionality of the i600 system.
- Any serious incident that has occurred related to the device should be reported to the manufacturer and the competent authority of the member state in which the user and patient are established.

5.2 Proper Training

WARNING

Before using your i600 system on patients:

- You should have been trained to use the system, or you should have read and fully understood this user guide.
 - You should be familiar with the safe use of the i600 system as detailed in this user guide.
 - Before use or after changing any settings, the user should check that the live image is displayed properly in the camera preview window of the program.
-

5.3 In Case of Equipment Failure

WARNING

If your i600 system is not working properly, or if you suspect that there is a problem with the equipment:

- Remove the device from the patient's mouth and discontinue use immediately.
- Disconnect the device from the PC and check for errors.
- Contact the manufacturer or authorized resellers.
- Modifications to the i600 system are prohibited by law as they may compromise the safety of the user, patient, or a third party.

5.4 Hygiene

WARNING

For clean working conditions and patient safety, ALWAYS wear clean surgical gloves when:

- Handling and replacing the tip.
- Using the i600 scanner on patients.
- Touching the i600 system.

The i600's main unit and its optical window should be kept clean at all times. Before using the i600 scanner on a patient, be sure to:

- Disinfect the i600 system
- Use a sterilized tip

5.5 Electrical Safety

WARNING

- The i600 system is a Class I device.
- To prevent electric shock, the i600 system must only be connected to a power source with a protective earth connection. If you are unable to insert the i600-supplied plug into the main outlet, contact a qualified electrician to replace the plug or outlet. Do not try to circumvent these safety guidelines.
- The i600 system only uses RF energy internally. The amount of RF radiation is low and does not interfere with surrounding electromagnetic radiation.
- There is a risk of electric shock if you attempt to access the inside of the i600 system. Only qualified service personnel should access the system.
- Do not connect the i600 system to a regular power strip or extension cord as these connections are not as safe as grounded outlets. Failure to adhere to these safety guidelines may result in the following hazards:
 - : The total short circuit current of all connected equipment may exceed the limit specified in EN / IEC 60601-1.
 - : The impedance of the ground connection may exceed the limit specified in EN / IEC 60601-1.
- Do not place liquids such as beverages near the i600 system and avoid spilling any liquid on the system.
- Condensation due to changes in temperature or humidity can cause moisture buildup inside the i600 unit, which may damage the system. Before connecting the i600 system to a power supply,

be sure to keep the i600 device at room temperature for at least two hours to prevent condensation. If condensation is visible on the product surface, the i600 should be left at room temperature for more than 8 hours.

- You should only disconnect the i600 system from the power supply via its power cord.
- When disconnecting the power cord, hold the plug surface to remove it.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services.
- Before disconnecting the power cord, make sure to turn off the power on the device using the power switch on the main unit.
- Only use the power adaptor supplied together with the i600. The use of other power adaptors may result in damage to the system.
- Avoid pulling on the communication cables, power cables, etc. used with the i600 system.

5.6 Eye Safety

WARNING

- The i600 system projects a bright light from its tip during scanning.
 - The bright light projected from the tip of the i600 is not harmful to the eyes. However, you should not look directly at the bright
-

light nor aim the light beam into the eyes of others. Generally, intense light sources can cause eyes to become brittle, and the likelihood of secondary exposure is high. As with other intense light source exposure, you may experience a temporary reduction in visual acuity, pain, discomfort, or visual impairment, all of which increase the risk of secondary accidents.

- Disclaimer for risks involving patients with epilepsy.
- The Medit i600 should not be used on patients that have been diagnosed with epilepsy due to the risk of seizures and injury. For the same reason, dental staff who have been diagnosed with epilepsy should not operate the Medit i600.

5.7 Explosion Hazards

WARNING

- The i600 system is not designed to be used near flammable liquids or gases, or in environments with high oxygen concentrations.
- There is a risk of explosion if you use the i600 system near flammable anesthetics.

5.8 Pacemaker and ICD Interference Risk

WARNING

- Do not use the i600 system on patients with pacemakers and ICD devices.
 - Check each manufacturer's instructions for interference by peripheral devices, such as computers used with the i600 system.
-

6 Electro-Magnetic Compatibility Information

6.1 Electro-Magnetic Emissions

This i600 is intended for use in the electromagnetic environment as specified below. The customer or the user of the i600 should ensure that it is used in such an environment.

| Emission Test | Compliance | Electromagnetic Environment - Guidance |
|--|------------|--|
| RF Emissions CISPR 11 | Group 1 | The i600 uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF Emissions CISPR 11 | Class A | The EUT is suitable for use in all establishments. This includes domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ Flicker emissions | Complies | |

Warning : This i600 is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the i600 or shielding the location.

6.2 Electro-Magnetic Immunity

Guidance 1

This i600 is intended for use in the electromagnetic environment as specified below. The customer or the user of the i600 system should ensure that it is used in such an environment.

| Immunity Test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|---|---|---|---|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with a synthetic material, a relative humidity of at least 30% is recommended. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV (for power supply lines) ± 1 kV (for input/output lines) | ± 2 kV (for power supply lines) ± 1 kV (for input/output lines) | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 0.5 kV, ± 1 kV differential mode ± 0.5 kV, ± 1 kV, ± 2 kV common mode | ± 0.5 kV, ± 1 kV differential mode ± 0.5 kV, ± 1 kV, ± 2 kV common mode | Mains power quality should be that of a typical commercial or hospital environment. |

| | | | |
|--|---|---|--|
| Voltage dips IEC 61000-4-11 | 0% U _T (100% dip in U _T) for 0.5 cycle at 50 Hz or 1 cycle at 60 Hz | 0% U _T (100% dip in U _T) for 0.5 cycle at 50 Hz or 1 cycle at 60 Hz | Mains power quality should be that of a typical commercial or hospital environment. If the user of the i600 image intensifier requires continued operation during power mains interruptions, it is recommended that the i600 image intensifier be powered from an uninterruptible power supply or a battery. |
| Short interruptions 61000-4-11 | 70% U _T (30% dip in U _T) for 20 cycles at 50 Hz or 30 cycle at 60 Hz | 70% U _T (30% dip in U _T) for 20 cycles at 50 Hz or 30 cycle at 60 Hz | |
| Voltage variations on power supply input lines 61000-4-11 | 0% U _T (100% dip in U _T) for 250 cycles at 50 Hz or 300 cycle at 60 Hz | 0% U _T (100% dip in U _T) for 250 cycles at 50 Hz or 300 cycle at 60 Hz | |
| Power frequency magnetic fields (50/60Hz) IEC 61000-4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a location in a typical commercial or hospital environment. |

NOTE : U_T is the main voltage (AC) prior to the application of the test level.

▪ Guidance 2

Recommended separation distances between portable and mobile communication equipment and the i600. The i600 is

intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the i600 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the i600 as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of the transmitter [W] | Separation distance according to the frequency of transmitter [m] | | | | |
|---|---|---|---|--|--|
| | IEC 60601 - 1 - 2: 2007 | | | IEC 60601 - 1 - 2: 2014 | |
| | 150 kHz to 80 MHz $d = 1.2\sqrt{P}$ | 800 MHz to 800 MHz $d = 1.2\sqrt{P}$ | 800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$ | 150 kHz to 80 MHz $d = 1.2\sqrt{P}$ | 80 MHz to 2.7 GHz $d = 2.0\sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 | 0.12 | 0.20 |
| 0.1 | 0.38 | 0.38 | 0.73 | 0.38 | 0.63 |
| 1 | 1.2 | 1.2 | 2.3 | 1.2 | 2.0 |
| 10 | 3.8 | 3.8 | 7.3 | 3.8 | 6.3 |
| 100 | 12 | 12 | 23 | 12 | 20 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

▪ **Guidance 3**

The i600 is intended for use in the electromagnetic environment specified below. The customer or the user of the i600 should assure that it is used in such an environment.

| Immunity Test | IEC 60601 test level | Compliance level | Recommended separation distance(d) | Electromagnetic environment - guidance |
|---------------|----------------------|------------------|------------------------------------|--|
|---------------|----------------------|------------------|------------------------------------|--|

| | | | | |
|----------------------------|--|-------|-------------------|---|
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz – 80 MHz Outside ISM Bands ^d 6 Vrms 150 kHz – 80 MHz In ISM Bands ^d | 3Vrms | $d = 1.2\sqrt{P}$ | Portable and mobile RF communications equipment, including cables, should be used no closer to any part of the i600 than the recommended separation distance as calculated using the equation below, according to the frequency of the transmitter. |
|----------------------------|--|-------|-------------------|---|

| | | |
|---------------------------|-------------------------|-------|
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2.7 GHz | 6 V/m |
|---------------------------|-------------------------|-------|

IEC 60601 - 1 - 2:2007 $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5 GHz

IEC 60601 - 1 - 2:2014 $d = 2.0\sqrt{P}$ 80 MHz to 2.7 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^{a)} should be less than the compliance level in each frequency range^{b)}
Interference may occur in the vicinity of equipment marked with the following symbol:

- NOTE 1 : At 80 MHz and 800 MHz, the higher frequency range applies.
- NOTE 2 : These guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.
 - a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the i600 is used exceeds the applicable RF compliance level above, the i600 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the i600.
 - b) When the frequency range exceeds 150 kHz – 80 MHz, the electric field strength should be not higher than 3 V/m.
 - c) The ISM (Industrial, Scientific, and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

▪ Guidance 4

The i600 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30cm (12 inches) to any part of the i600. Otherwise, degradation of the performance of this equipment could result.

| Immunity test | Band | Service | Modulation | IEC60601 test level | Compliance level |
|---|-----------------|---|--------------------------------|---------------------|------------------|
| Proximity fields From RF wireless Communications IEC61000 - 4 - 3 | 380 - 390 MHz | TETRA 400 | Pulse modulation 18 Hz | 27 V/m | 27 V/m |
| | 430 - 470 MHz | GMRS 460 FRS 460 | FM ±5 kHz deviation 1 kHz sine | 28 V/m | 28 V/m |
| | 704 - 787 MHz | LTE Band 13, 17 | Pulse modulation 217 Hz | 9 V/m | 9 V/m |
| | 800 - 960 MHz | GSM800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5 | Pulse modulation 18 Hz | 28 V/m | 28 V/m |
| | 1700 - 1990 MHz | GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1,2,4,25 UMTS | Pulse modulation 217 Hz | 28 V/m | 28 V/m |
| | 2400 - 2570 MHz | Bluetooth WLAN 802.11b/g/n RFID 2450 LTE Band 7 | Pulse modulation 217 Hz | 28 V/m | 28 V/m |
| | 5100 - 5800 MHz | WLAN 802.11a/n | Pulse modulation 217 Hz | 9 V/m | 9 V/m |

NOTE : If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m. The 1m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
 - b) The carrier shall be modulated using a 50% duty cycle square wave signal.
 - c) As an alternative to FM modulation, a 50% pulse modulation at 18 Hz may be used. While it does not represent actual modulation, it would be the worst case.
-

7 Specifications

| Model Name | MD-IS0100 |
|-----------------------------------|--|
| Trade Name | i600 |
| Rating | 9V $\overline{=}$, 3A |
| Applied part | Type BF |
| DC Adapter | |
| Model name | ATM036T-P120 |
| Input voltage | Universal 100-240 Vac / 50-60 Hz input, without any slide switch |
| Output | 12V $\overline{=}$, 3A |
| Case dimension | 100 x 50 x 33 mm (W x L x H) |
| EMI | CE / FCC Class B, Conduction & Radiation met |
| Protection | OVP (Over Voltage Protection) |
| | SCP (Short Circuit Protection) |
| | OCP (Over Current Protection) |
| Protection against electric shock | Class I |
| Mode of operation | Continuous |
| Handpiece | |
| Dimension | 248.2 x 44 x 47.4mm (W x L x H) |
| Weight | 231g |
| Power Hub | |
| Dimension | 68.2 x 31 x 14.9 mm (W x L x H) |
| Weight | 19 g |

| Calibration Tool | |
|-------------------------------------|---|
| Dimension | 124 x 54mm (H x Ø) |
| Weight | 220 g |
| Operating & Storage conditions | |
| Operating conditions | Temperature 18°C to 28°C |
| | Humidity 20 to 75% relative humidity (non-condensing) |
| | Air pressure 800 hPa to 1100 hPa |
| Storage conditions | Temperature -5°C to 45°C |
| | Humidity 20 to 80% relative humidity (non-condensing) |
| | Air pressure 800 hPa to 1100 hPa |
| Transport conditions | Temperature -5°C to 45°C |
| | Humidity 20 to 80% relative humidity (non-condensing) |
| | Air pressure 620 hPa to 1200 hPa |
| Emission limits per environment | |
| Environment | Hospital environment |
| Conducted and radiated RF EMISSIONS | CISPR 11 |

| | |
|----------------------------------|-------------------|
| Harmonic distortion | See IEC 61000-3-2 |
| Voltage Fluctuations and flicker | See IEC 61000-3-3 |

EC REP

EU representative

Meditrial Srl

Via Po 9 00198, Rome Italy

Email: ecrep@meditrial.eu Tel: +39-(0)6-45429780

Manufacturer

Address 1:23, Goryeodae-ro 22-gil, Seongbuk-gu, Seoul, 02855, Rep. of Korea

Address 2:F9, F10, 8, Yangpyeong-ro 25-gil, Yeongdeungpo-gu, Seoul, 07207, Rep. of Korea

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